

JORDAN TOF PRODUCTS, INC.

Formerly
R. M. Jordan Co., Inc.

TIME OF FLIGHT COMPONENTS ELECTRON GUN AND HIGH VOLTAGE PULSER PRODUCTS

[Printer Friendly Version
Of This Page](#)

[Reflectron TOF \(RETOF\)](#)

[Orthogonal Acceleration TOF
\(oaTOF\)](#)

[Ion Trap TOF](#)

[Linear Time Of Flight](#)

[Electron Gun Ionizers](#)

[Microchannel Plate Detectors](#)

[Home](#)

(Click image to enlarge)



Egun Mounted on 2.75 Conflat Flange

C-950 ELECTRON GUN IONIZER

The Electron Gun with its power supply can be used to inject electrons into a vacuum chamber. These electrons are primarily used to generate ions by electron impact. When used with the Pulser, it can be used as an EI source for a TOF Mass Spectrometer.

GENERAL DESCRIPTION

Electrons are emitted from a hot wire filament and accelerated to an energy level which is adjustable to between 40 and 100V. Most of these electrons are collected by the first lens element and are displayed on the meter as Emission Current. About 10% of them will pass through the slit and past the second lens element which is connected to the Focus Voltage. They then pass through the split lenses which are mechanically touching the Grid Plates of the Ion Source.

The Egun and Power Supply were designed to operate at the potential of the Grid Plates which can be as high as 3500VDC continuous. The Egun is mounted on a standard 2.75 CONFLAT flange.

The Filament Assembly is one inch in diameter and will extend into the vacuum as far as requested in order to make proper contact with the grid plates.

D-903 ELECTRON GUN POWER SUPPLY

EGUN POWER SUPPLY SPECIFICATIONS

EMISSION: 0 to 10 Milliamps

ELECTRON ENERGY: -40 to -100 Volts

FOCUS: 0 to -100 Volts

INPUT BIAS VOLTAGE: 3500 volts Maximum

CABINET SIZE 19.0" Width X 3.5" Height X 13.5" Depth

CABINET WEIGHT 9.5 Lbs.

INPUT POWER 100/120/220/240 Volts, 1 Phase, 50-60Hz

[C-950 Electron Gun and D-903 Power Supply Manual \(PDF\)](#)

(Click image to enlarge)



Electron Gun Filament Assembly

(Click image to enlarge)



D-903 Electron Gun Power Supply

(Click image to enlarge)



D-1040 High Voltage Pulser and D-1003 Power Supply

D-1040 HIGH VOLTAGE PULSER D-1003 PULSER POWER SUPPLY

In the usual TOF instrument, ions are generated with a Pulsed Laser which gives them all the same start time. In some applications ions are generated in a different way or a different place and there is a need for some way to inject them into the flight tube at the same time.

GENERAL DESCRIPTION

The Pulser is used to deliver a voltage pulse which can be varied in magnitude, duration and timing interval. These pulses when applied to elements of a TOF instrument can be used to extract bursts of ions and to gate or deflect Ion or Electron Beams.

It provides voltage steps of 11 to 400V from an internal power supply with a rise and fall time of <10 nanoseconds.

It can be biased up to 3500V continuous. This is adequate for an instrument with pulsed extraction voltage of less than 400V and acceleration voltage less than 3500V. Maximum repetition rate is 200 kHz.

The power supply can be wired for either positive or negative going pulses. It will also work with a positive or a negative bias voltage which is monitored on the front panel. NOTE: two wires must be changed to enable the meter to read upscale with a negative bias voltage. The power supply has a pulse delay circuit which

allows the user to adjust both delay and pulse duration. It must be triggered from an external pulse generator.

[D-1040 High Voltage Pulser and D-1003 Power Supply Manual \(PDF\)](#)

PULSER AND POWER SUPPLY SPECIFICATIONS

PULSE AMPLITUDE: +/-400 Volts from Input Bias Voltage
 BIAS VOLTAGE (NOT SUPPLIED): +/-3500V Maximum
 PULSE RISE/FALL TIME: <10nS
 PULSE DURATION: 50nS TO 12μS
 DELAY: 2μS to 1.2mS
 TRIGGER IN: +/- 1.3 to 10V TTL

CABINET SIZE 19.0" Width X 3.5" Height X 13.5" Depth
 CABINET WEIGHT 9.5 Lbs.
 REMOTE PULSER SIZE 5.5" Width X 1.1" Height X 4.3" Depth
 REMOTE PULSER WEIGHT 1 Lb.
 REMOTE PULSER OUTPUT SHV for direct mounting to flange feedthrough
 INPUT POWER 100/120/220/240 Volts, 1 Phase, 50-60Hz

[Click Images to Enlarge](#)



Description of Images and Product Information Downloads:

Electron Gun Ionizer for Linear and Reflectron Time of Flight Mass Spectrometers:

Electron Gun Mounted on 2.75 Inch Rotatable Conflat Flange (left)
 Electron Gun Mounted on C-870 Reflectron TOF Ion Source Assembly (right)

[EGUN Literature, Drawings and Wiring Diagram](#)
[EGUN Mount Detail for Shrouded Ion Source Assy](#)



Electron Gun Mounted on C-677 Linear Time of Flight Ion Source Assembly



Electron Gun Ionizer for Ion Trap Time of Flight Mass Spectrometers:

Electron Gun as Ionization Source for Ion Trap TOF (left)
 C-1550 Ion Trap and EGUN Assembly (right)

[Ion Trap EGUN Drawings and Wiring Diagram](#)



High Voltage Pulser and Power Supply:

D-1003 Pulser Power Supply (left)
 D-1040 High Voltage Pulser (right)

[Click Images to Enlarge](#)


[Electron Gun and Pulser Price List](#)

User Manuals:

[C-950 Electron Gun and D-903 Power Supply Manual \(PDF\)](#) Filesize: approx. 350 KB.

[D-1040 High Voltage Pulser and D-1003 Power Supply Manual \(PDF\)](#) Filesize: approx. 450 KB.

Click on link to open. To download, right click on link and "save target as".



[Home](#)

[Product Catalog](#)

[Time Of Flight Tutorial](#)

[Price List](#)

[Contact Us](#)

Jordan TOF Products, Inc.

990 Golden Gate Terrace
 Grass Valley, CA 95945
 530-272-4580 / 530-272-2955 [fax]
 E-mail: info@rmjordan.com

Copyright © Jordan TOF Products, Inc. All rights reserved.
 Thank you for visiting <http://www.rmjordan.com/egun-tof.html>

Electron Gun and High Voltage Pulser Time Of Flight / TOF / Mass Spectrometer Components